

Abbott Issues Voluntary Recall Of Certain FreeStyle® And FreeStyle Lite® Blood Glucose Test Strips In The United States

Abbott Park, Illinois — Abbott today announced it is initiating a voluntary recall of 20 lots of FreeStyle® and FreeStyle Lite® Blood Glucose Test Strips in the United States. These lots of test strips may produce erroneously low blood glucose results when used with both "FreeStyle® Blood Glucose Meter" and "FreeStyle® Flash Blood Glucose Meter" [neither of which have been in production since 2010], as well as the OmniPod® Insulin Management System¹.

When the test strips are used with the newer FreeStyle brand meters including FreeStyle® Freedom Blood Glucose Meter, FreeStyle Lite® Blood Glucose Meter* and FreeStyle Freedom® Lite Blood Glucose Meter*, the blood glucose test results are not affected. Testing with the FreeStyle InsuLinx® Blood Glucose Meter is not affected by this action, as FreeStyle InsuLinx Blood Glucose Meter uses FreeStyle InsuLinx test strips.

The affected test strips are from lots that hold expiration dates from between May 2014 and March 2015, and can be identified as follows:

Lot Numbers		
1281732	1363321	1367917
1283345	1365056	1373262
1283603	1365920	1374907
1285007	1365934	1366515
1366111	1366337	1366006
1363015	1366347	1350414
1363109	1365921	

The company is notifying healthcare professionals, pharmacies, distributors and customers about the recall; customers affected by this action are instructed to call Abbott's diabetes care customer service at 1-888-736-9869 for a replacement of the affected test strips at no charge.

For more information, please visit abbottdiabetescare.com.

About Abbott Diabetes Care

Abbott Diabetes Care, based in Alameda, Calif., is a leader in developing, manufacturing and marketing glucose monitoring systems designed to help people better manage their diabetes. Additional information about Abbott Diabetes Care may be found at www.abbottdiabetescare.com.

About Abbott

Abbott (NYSE: ABT) is a global healthcare company devoted to improving life through the development of products and technologies that span the breadth of healthcare. With a portfolio of leading, science-based offerings in diagnostics, medical devices, nutritionals and branded generic pharmaceuticals, Abbott serves people in more than 150 countries and employs approximately 70,000 people.

Visit Abbott at www.abbott.com and connect with us on Twitter at @AbbottNews.

¹ OmniPod is a registered trademark of Insulet Corporation

* FreeStyle Lite Blood Glucose Test Strips are intended to be used with FreeStyle Lite and FreeStyle Freedom Lite meters only

Media:

Jonathon Hamilton
Vicky Assardo
(408) 624-0314
(510) 207-7623

Financial:

Tina Ventura
(847) 935-9390
